

# Silensor-sl

Verbinderswechsel  
exchange of a connector



R<sub>x</sub> only









23°F  
-5°C



45°C  
113°F

• changement d'attache • cambio de conector • sostituzione congiuntore • substituição do conector




- loser Anker
- loose anchor
- ancre lâche
- ancla suelta
- ancoraggio staccato

• soltar suporte





- Anker halten
- hold anchor
- tenir ancre
- mantener ancla
- tenere premuto l'ancoraggio
- segurar suporte




restliche Verbinder hier aufkleben  
glue remaining connectors here  
coller attaches restantes ici  
pegar conectores restantes aquí  
colla premuto l'ancoraggio qui  
colar aqui os restantes conectores



zuerst verwendeter Verbinder  
first used connector  
attache en premier utilisée  
conector utilizó por primera vez  
congiuntore prima utilizzazione  
conector a utilizar na primeira vez



oberer Eckzahn  
upper canine  
canine supérieure  
canino superior  
canino superiore  
canino superior



unterer Molarenbereich  
lower molar area  
région molaire inférieure  
zona de molares inferior  
zona dei molari inferiori  
zona do molar inferior



sl-Zange/sl-pliers/sl-pince/sl-pinza  
Empfehlung zum leichteren Verbinderwechsel/  
Recommendation for an easier connector  
replacement, Ref.: 59 60 60

